

Whole Exome Sequencing

Gene package Ciliopathy, version 6, 21-2-2020



Technical information

DNA was enriched using Agilent SureSelect DNA + SureSelect OneSeq 300kb CNV Backbone + Human All Exon V7 capture and paired-end sequenced on the Illumina platform (outsourced). The aim is to obtain 10 Giga base pairs per exome with a mapped fraction of 0.99. The average coverage of the exome is ~50x. Duplicate and non-unique reads are excluded. Data are demultiplexed with bcl2fastq Conversion Software from Illumina. Reads are mapped to the genome using the BWA-MEM algorithm (reference: <http://bio-bwa.sourceforge.net/>). Variant detection is performed by the Genome Analysis Toolkit HaplotypeCaller (reference: <http://www.broadinstitute.org/gatk/>). The detected variants are filtered and annotated with Cartagenia software and classified with Alamut Visual. It is not excluded that pathogenic mutations are being missed using this technology. At this moment, there is not enough information about the sensitivity of this technique with respect to the detection of deletions and duplications of more than 5 nucleotides and of somatic mosaic mutations (all types of sequence changes).



Dept. Clinical Genetics

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
ACVR2B	602730	105	100	100	98
ADAMTS9	605421	67	99	96	90
AHI1	608894	67	100	98	90
ALMS1	606844	100	100	100	99
ANKS6	615370	89	95	93	90
ARL13B	608922	63	100	100	94
ARL6	608845	43	96	95	83
ARMC4	615408	67	98	96	92
ARMC9	617612	76	100	98	93
B9D1	614144	90	100	100	100
B9D2	611951	82	100	100	100
BBIP1	613605	45	100	96	86
BBS1	209901	126	100	100	100
BBS10	610148	66	100	100	100
BBS12	610683	59	100	100	99
BBS2	606151	77	100	100	97
BBS4	600374	76	100	100	95

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
BBS5	603650	65	100	98	89
BBS7	607590	64	100	99	95
BBS9	607968	59	96	95	90
C2CD3	615944	80	100	100	99
C8orf37	614477	78	100	100	96
CC2D2A	612013	66	100	100	95
CCDC103	614677	83	100	100	100
CCDC114	615038	125	100	100	100
CCDC151	615956	104	100	100	100
CCDC28B	610162	83	100	100	100
CCDC39	613798	80	100	99	94
CCDC40	613799	100	100	100	100
CCDC65	611088	57	100	100	94
CCNO	607752	113	100	100	100
CENPF	600236	70	100	99	95
CEP104	616690	71	100	98	93
CEP120	613446	75	100	100	96
CEP164	614848	93	100	100	98
CEP290	610142	71	100	98	89
CEP41	610523	69	100	100	94
CEP55	610000	71	100	100	97
CEP83	615847	60	100	96	81
CFAP298	615494	75	100	100	96
CFAP300	618058	53	100	98	87
CFAP410	603191	141	100	100	100
CFAP53	614759	98	100	100	98
CFC1	605194	70	41	34	28
CPLANE1	614571	71	100	99	95
CSPP1	611654	83	100	100	97
DCDC2	605755	96	100	100	97
DDX59	615464	62	100	98	92
DNAAF1	613190	131	100	100	99
DNAAF2	612517	124	100	100	100
DNAAF3	614566	98	100	100	99
DNAAF4	608706	54	100	99	89

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
DNAAF5	614864	104	96	86	79
DNAH1	603332	127	100	100	100
DNAH11	603339	74	100	100	96
DNAH17	610063	94	100	99	96
DNAH5	603335	74	100	100	96
DNAH9	603330	93	100	100	98
DNAI1	604366	91	100	98	97
DNAI2	605483	131	100	100	99
DNAJB13	610263	88	100	100	100
DNAL1	610062	77	100	100	95
DNHD1	617277	117	100	100	100
DRC1	615288	81	100	100	98
DYNC2H1	603297	71	100	99	91
DYNC2LI1	617083	53	100	100	91
EVC	604831	98	96	94	94
EVC2	607261	87	100	99	95
EXOC8	615283	118	100	100	100
FAM149B1	618413	61	100	93	83
FOXJ1	602291	116	100	100	100
FUZ	610622	109	100	100	100
GAS2L2	611398	154	100	100	100
GAS8	605178	85	99	97	93
GDF1	602880	69	100	93	81
GLIS2	608539	117	100	100	100
HYDIN	610812	46	82	69	57
HYLS1	610693	74	100	100	100
IFT122	606045	111	100	100	99
IFT140	614620	113	100	100	98
IFT172	607386	75	100	100	96
IFT27	615870	84	100	100	99
IFT43	614068	75	100	100	100
IFT52	617094	50	100	99	91
IFT80	611177	62	100	98	87
IFT81	605489	47	99	94	82
INPP5E	613037	113	100	100	97

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
INTU	610621	59	100	99	92
INVS	243305	81	100	100	97
IQCB1	609237	54	100	94	78
KCTD3	613272	61	100	100	96
KIAA0556	616650	97	100	99	97
KIAA0586	610178	72	100	98	94
KIAA0753	617112	52	100	98	91
KIF14	611279	74	100	99	91
KIF7	611254	107	98	96	93
LCA5	611408	74	100	100	99
LEFTY2	601877	144	100	100	99
LRRC56	618227	118	100	100	98
LRRC6	614930	79	100	99	93
LZTFL1	606568	79	100	100	93
MAPKBP1	616786	108	100	100	98
MCIDAS	614086	97	100	100	92
MKKS	604896	80	100	100	100
MKS1	609883	117	100	100	99
MMP21	608416	81	100	96	88
NCAPG2	608532	51	100	94	79
NEK1	604588	56	100	98	89
NEK8	609799	132	100	100	100
NME8	607421	65	100	100	93
NODAL	601265	119	100	100	100
NPHP1	607100	58	100	98	88
NPHP3	608002	66	100	99	94
NPHP4	607215	110	100	100	100
OCRL	300535	43	100	95	80
OFD1	300170	43	100	95	75
PDE6D	602676	86	100	100	100
PIBF1	607532	42	100	93	72
PIH1D3	300933	40	100	92	66
PIK3C2A	603601	53	100	98	91
PKD1	601313	93	97	95	90
PKD2	173910	62	100	97	89

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
PKHD1	606702	82	100	100	98
POC1A	614783	90	100	100	100
RPGRIP1	605446	88	100	100	96
RPGRIP1L	610937	60	98	96	89
RSPH1	609314	80	100	100	95
RSPH3	615876	89	100	99	94
RSPH4A	612647	100	100	100	100
RSPH9	612648	125	100	100	100
SCLT1	611399	51	100	97	87
SDCCAG8	613524	92	100	100	96
SPAG1	603395	64	100	99	91
SPATA7	609868	62	100	99	94
SUFU	607035	104	100	100	100
TBC1D32	615867	61	100	98	90
TCTEX1D2	617353	89	100	100	93
TCTN1	609863	94	100	100	98
TCTN2	613846	84	100	100	97
TCTN3	613847	66	100	100	97
TMEM107	616183	100	100	100	94
TMEM138	614459	54	100	100	100
TMEM216	613277	87	100	100	95
TMEM231	614949	96	100	95	90
TMEM237	614423	57	100	98	87
TMEM260	617449	58	100	98	93
TMEM67	609884	77	100	100	94
TRAF3IP1	607380	55	100	97	86
TRIM32	602290	107	100	100	100
TTBK2	611695	74	100	100	99
TTC21B	612014	75	100	99	94
TTC25	617095	71	100	99	97
TTC26	617453	55	100	99	93
TTC8	608132	66	100	100	93
TULP1	602280	108	100	100	100
WDPCP	613580	72	100	99	93
WDR19	608151	71	100	100	97

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	median depth	% covered >10x	% covered >20x	% covered >30x
WDR34	613363	116	100	100	100
WDR35	613602	62	100	100	94
WDR60	615462	80	100	99	94
XPNPEP3	613553	63	100	100	96
ZIC3	300265	101	100	100	97
ZMYND10	607070	130	100	100	100
ZNF423	604557	166	100	100	100

- OMIM release used: 8-9-2019

- The statistics above are based on a set of 100 samples

- Median depth is the median of the mean sequence depth over the protein coding exons (± 10 bp flanking introns) of the longest transcript

- % Covered 10x , 20x and 30 x describes the percentage of a gene's coding sequence (± 10 bp flanking introns) that is covered at least 10x, 20x or 30x