

# Whole Exome Sequencing

## Gene package Cardiomyopathy, version 10, 26-2-2021



### Technical information

DNA was enriched using Agilent SureSelect DNA + SureSelect OneSeq 300kb CNV Backbone + Human All Exon V7 capture and paired-end sequenced on the Illumina platform (outsourced). The aim is to obtain 10 Giga base pairs per exome with a mapped fraction of 0.99. The average coverage of the exome is ~50x. Duplicate and non-unique reads are excluded. Data are demultiplexed with bcl2fastq Conversion Software from Illumina. Reads are mapped to the genome using the BWA-MEM algorithm (reference: <http://bio-bwa.sourceforge.net/>). Sequence variant detection is performed by the Genome Analysis Toolkit HaplotypeCaller (reference: <http://www.broadinstitute.org/gatk/>). The detected sequence variants are filtered and annotated with Alissa Interpret software and classified with Alamut Visual. Copy variant detection is performed using the BAM multiscale reference method using depth of coverage analysis and dynamical bins in NexusClinical. The detected copy number variants are filtered and annotated with the NexusClinical software and classified using UCSC Genome Browser (NCBI37/hg19). It is not excluded that pathogenic variants are being missed using this technology. At this moment, there is not enough information about the sensitivity of this technique with respect to the detection of deletions and duplications of more than 5 nucleotides and of somatic mosaic mutations (all types of sequence changes).



**Dept. Clinical Genetics**

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	% covered $\geq 10x$	% covered $\geq 20x$	% covered $\geq 30x$	% covered $\geq 50x$
ACAD9	611103	100	100	100	98.36
ACTC1	102540	100	100	100	95.28
ACTN2	102573	100	100	99.61	94.81
ALPK3	617608	96.10	94.47	93.09	83.40
ANO5	608662	100	99.31	96.43	89.27
BAG3	603883	100	100	98.64	94.82
CASQ2	114251	100	100	99.93	97.28
CDH2	114020	98.80	97.37	95.26	89.38
CRYAB	123590	100	100	100	100
CSRP3	600824	100	100	100	100
DES	125660	100	100	99.34	91.76
DSC2	125645	100	97.12	96.59	93.52
DSG2	125671	99.44	97.52	96.21	91.68
DSP	125647	100	100	100	99.04
DTNA	601239	100	100	100	100
EMD	300384	100	96.66	92.01	76.59
EYA4	603550	100	100	100	98.36
FHL1	300163	100	100	98.83	94.26

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	% covered ≥10x	% covered ≥20x	% covered ≥30x	% covered ≥50x
FHOD3	609691	100	98.93	96.35	88.12
FKTN	607440	100	100	99.39	91.81
FLNC	102565	100	99.99	99.23	94.11
GATAD1	614518	100	100	98.90	91.73
GLA	300644	100	100	99.72	95.09
HCN4	605206	100	99.50	96.83	83.07
ILK	602366	99.94	99.94	99.94	99.94
JPH2	605267	100	100	100	100
JUP	173325	100	100	99.80	94.75
KIF20A	605664	100	100	100	99.60
LAMP2	309060	100	99.96	94.83	71.09
LMNA	150330	100	100	96.42	75.25
MIB1	608677	100	99.21	96.18	91.67
MYBPC3	600958	100	99.36	97.28	93.26
MYH6	160710	100	99.80	98.65	93.05
MYH7	160760	100	100	99.37	93.53
MYL2	160781	100	100	100	98.51
MYL3	160790	100	100	100	78.44
MYLK2	606566	100	100	100	96.13
MYOZ1	605603	100	100	100	96.44
MYOZ2	605602	100	100	99.33	93.83
MYPN	608517	100	100	99.49	95.01
NEBL	605491	100	100	100	99.21
NEXN	613121	100	94.24	84.90	68.32
NKX2-5	600584	100	97.26	86.00	31.85
PDLIM3	605889	100	100	100	99.59
PKP2	602861	100	99.37	94.77	90.56
PLN	172405	100	100	100	98.86
PPA2	609988	100	93.63	82.13	64.78
PPCS	609853	100	100	97.63	85.20
PPP1R13L	607463	100	94.22	81.13	63.12
PRDM16	605557	100	98.39	96.28	90.31
PRKAG2	602743	97.51	97.11	92.70	84.41
PTPN11	176876	98.37	98.37	98.37	98.28
RAF1	164760	100	100	97.97	93.07
RBM20	613171	100	100	99.18	90.03
RYR2	180902	100	99.37	97.91	94.91
SCN5A	600163	100	99.89	99.07	95.89
SGCD	601411	100	100	100	92.68
SOS1	182530	100	97.71	97.35	91.97

HGNC approved gene symbol	OMIM gene ID (active link to omim.org)	% covered ≥10x	% covered ≥20x	% covered ≥30x	% covered ≥50x
TAZ	300394	100	100	98.38	90.52
TBX20	606061	100	100	100	94.64
TCAP	604488	100	100	100	90.94
TGFB3	190230	100	100	100	96.29
TMEM43	612048	100	100	99.93	96.32
TNNC1	191040	100	100	100	98.01
TNNI3	191044	100	100	100	87.09
TNNI3K	613932	100	100	100	97.39
TNNT2	191045	100	100	98.81	87.81
TPM1	191010	99.84	91.83	83.55	81.11
TRIM63	606131	100	100	100	96.53
TTN	188840	100	100	99.94	99.55
TTR	176300	100	100	100	100
TXNRD2	606448	99.69	97.31	93.57	86.82
VCL	193065	100	99.93	98.50	91.93

- OMIM release used: 18-2-2021

- The statistics above are based on a set of 100 samples

- % Covered 10x , 20x, 30x and 50x describes the percentage of a gene's coding sequence (±10bp flanking introns) that is covered at least 10x, 20x, 30x or 50x