

Whole Exome Sequencing

Gene package Very early onset inflammatory bowel disease / congenital diarrhea, version 5, 26-2-2021



Technical information

DNA was enriched using Agilent SureSelect DNA + SureSelect OneSeq 300kb CNV Backbone + Human All Exon V7 capture and paired-end sequenced on the Illumina platform (outsourced). The aim is to obtain 10 Giga base pairs per exome with a mapped fraction of 0.99. The average coverage of the exome is ~50x. Duplicate and non-unique reads are excluded. Data are demultiplexed with bcl2fastq Conversion Software from Illumina. Reads are mapped to the genome using the BWA-MEM algorithm (reference: <http://bio-bwa.sourceforge.net/>). Sequence variant detection is performed by the Genome Analysis Toolkit HaplotypeCaller (reference: <http://www.broadinstitute.org/gatk/>). The detected sequence variants are filtered and annotated with Alissa Interpret software and classified with Alamut Visual. Copy variant detection is performed using the BAM multiscale reference method using depth of coverage analysis and dynamical bins in NexusClinical. The detected copy number variants are filtered and annotated with the NexusClinical software and classified using UCSC Genome Browser (NCBI37/hg19). It is not excluded that pathogenic variants are being missed using this technology. At this moment, there is not enough information about the sensitivity of this technique with respect to the detection of deletions and duplications of more than 5 nucleotides and of somatic mosaic mutations (all types of sequence changes).



Dept. Clinical Genetics

HGNC approved gene symbol	OMIM gene ID	% covered $\geq 10x$	% covered $\geq 20x$	% covered $\geq 30x$	% covered $\geq 50x$
ACTG2	102545	100	100	100	99.30
ADAM17	603639	100	98.56	95.42	88.36
AIRE	607358	100	99.35	95.72	87.42
ALPI	171740	100	100	100	97.39
ANKZF1	617541	100	99.75	96.35	84.80
ARPC1B	604223	100	100	100	95.56
ARX	300382	85.61	80.04	71.89	58.82
BACH2	605394	100	100	100	98.21
CARD8	609051	100	100	100	96.43
CARMIL2	610859	97.81	94.63	92.52	82.62
CD19	107265	100	99.56	97.49	90.21
CD40	109535	100	100	100	98.05
CD40LG	300386	100	100	100	96.60
CD55	125240	100	96.43	91.53	78.51
CLMP	611693	100	100	100	96.19
COG6	606977	100	100	99.53	92.92
CR2	120650	100	100	99.58	95.80
CTLA4	123890	100	100	100	99.33
CYBA	608508	94.68	83.76	65.53	50.71
CYBB	300481	100	99.95	97.03	83.25
CYBC1	No ID	100	96.84	89.21	83.55
DEF6	610094	100	96.55	92.66	87.20
DGAT1	604900	95.59	90.52	87.80	86.72
DOCK2	603122	100	99.54	98.91	95.43

HGNC approved gene symbol	OMIM gene ID	% covered ≥10x	% covered ≥20x	% covered ≥30x	% covered ≥50x
DOCK8	611432	100	99.53	98.26	93.32
DUOX2	606759	100	97.09	95.22	87.27
EGFR	131550	100	99.77	98.29	94.11
ELF4	300775	100	99.86	98.53	92.18
EPCAM	185535	100	98.75	92.47	90.73
FOXP3	300292	100	100	99.17	89.13
GPIHBP1	612757	100	100	100	91.69
GUCY2C	601330	100	99.79	99.08	93.67
HPS1	604982	100	100	99.63	88.78
HPS4	606682	100	100	98.64	92.29
HPS6	607522	100	100	95.65	78.12
ICOS	604558	100	100	100	100
IKBKG	300248	33.32	25.95	25.95	22.46
IL10	124092	100	100	100	99.13
IL10RA	146933	100	100	98.56	94.08
IL10RB	123889	100	94.66	93.72	89.42
IL21	605384	100	99.15	92.92	83.28
IL21R	605383	100	100	100	95.94
IL2RA	147730	100	100	100	98.51
IL2RB	146710	100	100	98.31	90.97
ITCH	606409	98.02	96.02	93.74	88.49
ITGB2	600065	100	98.83	96.43	85.40
LCT	603202	100	100	99.55	98.08
LRBA	606453	99.80	99.46	98.51	95.25
MALT1	604860	95.57	91.25	90.18	88.28
MEFV	608107	100	100	99.10	92.94
MPI	154550	100	100	100	99.16
MVK	251170	100	100	100	96.47
MYO5B	606540	100	99.86	98.93	94.24
NCF1	608512	58.99	51.83	39.42	24.06
NCF2	608515	100	100	100	98.11
NCF4	601488	100	100	98.40	83.32
NEUROG3	604882	100	100	96.68	90.94
NFKBIA	164008	100	100	100	93.18
NLRC4	606831	100	100	100	99.81
NOD2	605956	100	100	99.78	96.90
OTULIN	615712	91.89	85.62	85.62	80.64
PCSK1	162150	100	100	99.67	96.49
PLA2G4A	600522	100	100	99.44	96.60
PLCG2	600220	100	100	99.33	93.91
PLVAP	607647	100	100	98.24	89.87
PNLIP	246600	100	100	99.88	96.42
POLA1	312040	100	98.83	92.93	75.88
PTPN2	176887	98.52	93.86	93.38	87.03

HGNC approved gene symbol	OMIM gene ID	% covered ≥10x	% covered ≥20x	% covered ≥30x	% covered ≥50x
RAC2	602049	100	100	100	98.42
RFX6	612659	100	100	99.29	96.29
RIPK1	603453	100	100	98.60	92.41
SAR1B	607690	99.30	90.13	87.96	87.96
SH2D1A	300490	100	100	100	89.33
SI	609845	99.97	98.07	95.42	89.62
SKIV2L	600478	100	100	98.56	92.86
SLC10A2	601295	100	100	100	94.54
SLC26A3	126650	100	100	100	96.21
SLC2A2	138160	100	100	100	96.07
SLC37A4	602671	99.87	99.87	99.10	93.11
SLC39A4	607059	100	100	100	96.45
SLC5A1	182380	100	100	98.76	93.26
SLC7A7	603593	100	100	100	97.31
SLC9A3	182307	100	98.06	94.88	84.38
SPINT2	605124	100	100	95.26	81.03
STAT1	600555	97.05	97.05	97.05	96.66
STAT3	102582	100	100	99.75	94.86
STAT5B	604260	100	98.35	96.03	87.06
STX3	600876	100	100	95.31	91.19
STXBP2	601717	100	97.15	92.77	77.46
TCN2	613441	100	100	99.90	95.24
TGFB1	190180	100	100	100	99.01
TMPRSS15	606635	96.29	88.28	87.84	84.97
TNFAIP3	191163	100	100	99.01	96.82
TNFRSF13B	604907	100	100	98.77	90.35
TNFRSF13C	606269	91.34	66.50	60.13	51.23
TOM1	604700	100	100	100	97.16
TRIM22	606559	100	100	98.23	87.85
TTC37	614589	100	100	99.48	97.36
TTC7A	609332	100	100	99.05	90.00
TYMP	131222	100	100	94.91	78.46
WAS	300392	100	100	98.42	83.59
WIPF1	602357	100	98.97	96.18	77.02
WNT2B	601968	100	100	99.90	96.83
XIAP	300079	100	100	95.53	84.45
ZAP70	176947	100	99.67	97.76	91.34

- OMIM release used: 18-2-2021

- The statistics above are based on a set of 100 samples

- % Covered 10x , 20x, 30x and 50x describes the percentage of a gene's coding sequence (±10bp flanking introns) that is covered at least 10x, 20x, 30x or 50x